
SPECIFICATION SHEET

| | | | |
|-----------------------|--|-------------------|----------|
| Product Name: | McAb-HIV-1P24 | Catalog #: | L2H00503 |
| Description: | MAb to HIV-1 p24 Monoclonal Antibody to Human Immunodeficiency Virus Type 1 (HIV-1), p24 specific | | |
| Specificity: | HIV-1 p24 (Recombinant antigen) | | |
| Clone: | MLP3E12 | | |
| Host Animal: | Mouse | | |
| Isotype: | IgG | | |
| Source: | Ascites | | |
| Immunogen: | Recombinant Antigen | | |
| Format: | Purified, Liquid | | |
| Concentration: | 1.8mg/ml (OD280nm) | Lot: | L190401 |
| Buffer: | 0.01M PBS, pH 7.4 | | |
| Preservative: | None | | |
| Applications: | Suitable as a detection antibody in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. | | |
| Storage: | Store at -20°C. Avoid multiple freeze/thaw cycles. | | |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.