

## SPECIFICATION SHEET

<b>Product Name:</b>	Mouse anti human HAMA	<b>Catalog #:</b> L1I006
<b>Description:</b>	Mouse Antibody to human HAMA	
<b>Specificity:</b>	human HAMA	
<b>Host Animal:</b>	Mouse	
<b>Immunogen:</b>	Purified HAMA	
<b>Format:</b>	Liquid	
<b>Purification:</b>	>80% pure (SDS-PAGE).	
<b>Concentration:</b>	45mg/ml	<b>Lot.:</b> L190801
<b>Buffer:</b>	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2	
<b>Preservative:</b>	None	
<b>Applications:</b>	Specific methodologies have not been tested using this product.	
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.	

**FOR RESEARCH USE ONLY**