



Hainan Medikamende Trading Co., Ltd.

海南省迈迪卡门德商贸有限责任公司  
Hainan Medikamende Trading Co., Ltd.

---

## SPECIFICATION SHEET

**Description:** Recombinant human HE4 Antigen **Catalog #:** L2C02601

**Molecular Weight:** 14kd (calculated from sequence)  
14KD(determined by SDS gelelectrophoresis)

**Source:** E. coli

**Format:** Purified, Liquid

**Purification:** > 90% (determined by SDS gelelectrophoresis)

**Concentration:** 1mg/ml (OD280 OD260) **Lot:** L160701

**Buffer:** 20Mm Tris-HCL, pH 8.0 + 0.2

**Preservative:** 0.02%NaN3

**Storage:** Store at -20 C for 3years. Avoid multiple freeze/thaw cycles.

---

海南省迈迪卡门德商贸有限责任公司  
地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502  
邮编：572000  
电话：+86-898-31320087，传真：+86-898-31320087  
网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.  
Adress: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,  
Xiuyingqu, Haikou, 572000, Hainan, P.R. China  
Phone: +86-898-31320087 Fax: +86-898-31320087  
Website: www.medical-ivd.com Emai: info@medical-ivd.com