



SPECIFICATION SHEET

Catalog #:	L2C01601	Lot #:	L180301
Description:	Recombinant human β 2-MG Antigen (B2M) Recombinant human Beta 2 Microglobulin Antigen		
Molecular Weight:	17kd (calculated from sequence) 17KD(determined by SDS gelelectrophoresis)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	> 90% (determined by SDS gelelectrophoresis)		
Concentration:	2.91mg/ml		
Buffer:	20mM Tris, pH 8.0 \pm 0.2, 1mm EDTA		
Preservative:	None		
Storage:	Store at -20 C. Avoid multiple freeze/thaw cycles.		