

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Description: L-T₃ 3,3', 5-Triiodo-L-thyronine Sodium Salt

Catalog #: L2G01001

Molecular Formula: C₁₅H₁₁I₃NNaO₄

MW: 672.96

Purity: 95%

Lot: L160601

Forma: Powder

Storage: Store at 2 ~ 8°C, Against light or wet.

Dissolved: Dissolve one by one drops of NaOH 0.1mol/L, and make it dissolve completely. The solution is diluted to the appropriate concentration with 50% propylene glycol solution.(Note: this method is only applicable to the calibration of the product, for reference. Other uses of the dissolution method please refer to the relevant information.)