

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>ProductName:</b>	McAb to S100-A8	<b>Catalog #:</b>	L6B00401
<b>Description:</b>	Mouset anti Recombinant Human S100 Calcium Binding Protein A8		
<b>Specificity:</b>	Recognizes an epitope on Human S100 Calcium Binding Protein A8		
<b>Immunogen:</b>	Recombinant Human S100 Calcium Binding Protein A8		
<b>Host Animal:</b>	Mouse		
<b>Source:</b>	Serum		
<b>Format:</b>	Purified, Liquid		
<b>Concentration:</b>	4.9mg/ml	<b>LOT:</b>	L180801
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Storage&amp;Shelflife:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		