

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	McAb to AMH	Catalog #:	L1F00801				
Description:	Monoclonal Antibody to Human anti-Mullerian hormone, AMH						
Clone:	AM001						
Host Animal:	Mouse						
Isotype:	IgG						
Source:	Ascites						
Immunogen:	Recombinant Human anti-Mullerian hormone, AMH						
Format:	Purified, Liquid						
Concentration&Lot:	4.5 mg/ml (OD280nm)	Lot #:	L180401				
Buffer:	10mM PBS pH 7.4 ± 0.2						
Preservative:	None						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1F00801</td><td>L1F00802, L1F00803</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1F00801	L1F00802, L1F00803
<u>Capture</u>	<u>Detection</u>						
L1F00801	L1F00802, L1F00803						
Storage:	Store at -20°C for at least 3 years.						