

SPECIFICATION SHEET

Product Name:	McAb to Cancer Antigen CA125	Catalog #:	L1C01202				
Description:	Monoclonal Antibody to Cancer Antigen CA125 (CA125), Ovarian Cancer						
Specificity:	Recognizes an epitope on the CA125 containing mucin. The antibody does recognize the M11-like epitope (Group B1), And pairs well with the OC-125 specific monoclonal antibodies when used in a sandwich assay.						
Clone:	CA002						
Host Animal:	Mouse						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Natural Human CA125 Antigen						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration&Lot:	4mg/ml	Lot :	L20200401				
Buffer:	0.15M NaCl, 10mM PB,PH7.5						
Preservative:	0.05%NaN ₃						
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C01201</td><td>L1C01202</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C01201	L1C01202
<u>Capture</u>	<u>Detection</u>						
L1C01201	L1C01202						
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY.