



Hainan Medikamende Trading Co., Ltd.

海南省迈迪卡门德商贸有限责任公司
Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Product Name: McAb to Human Ferritin **Catalog #:** L1C01502

Description: Monoclonal Antibody to Human Ferritin

Specificity: Human Ferritin

Clone: FE002

Host Animal: Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice.

Source: Ascites

Immunogen: Purified human ferritin

Format: Purified, Liquid

Purification: Protein G chromatography

Concentration&Lot: 5.4mg/ml **Lot :** L180709

Buffer: 0.01M PBS, 0.15M NaCl, pH 7.4

Preservative: None

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Recommended pair for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
L1C01501	L1C01502

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.

海南省迈迪卡门德商贸有限责任公司
地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
邮编：572000
电话：+86-898-31320087，传真：+86-898-31320087
网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.
Address: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,
Xiuyingqu, Haikou, 572000, Hainan, P.R. China
Phone: +86-898-31320087 Fax: +86-898-31320087
Website: www.medical-ivd.com Emai: info@medical-ivd.com