

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00801	Lot #: L180801				
Description:	MAb to FSH Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH) alpha subunit					
Clone:	FS001					
Host Animal:	Mouse					
Isotype:	IgG ₁ ,					
Source:	Ascites					
Immunogen:	Purified natural FSH alpha					
Format:	Purified, Liquid					
Concentration:	1.7 mg/ml					
Buffer:	10mM PBS pH 7.4 ± 0.2					
Preservative:	None					
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00801</td><td>L1G00802</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00801	L1G00802
<u>Capture</u>	<u>Detection</u>					
L1G00801	L1G00802					
Storage:	Short term (up to 2 Weeks) store at 2–8°C. Long term, store at <-20°C.					