

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

| | | | | | | |
|-----------------------|--|-----------------------|----------------|------------------|----------|----------|
| Catalog #: | L1G00102 | Lot #: L160701 | | | | |
| Description: | McAb to PRL Monoclonal Antibody to prolactin, PRL. | | | | | |
| Specificity: | Human prolactin, PRL | | | | | |
| Clone: | PR002 | | | | | |
| Host Animal: | Mouse | | | | | |
| Source: | Ascites | | | | | |
| Immunogen: | Purified Recombinant prolactin | | | | | |
| Format: | Purified, Liquid | | | | | |
| Purification: | ≥90% pure (SDS-PAGE). | | | | | |
| Concentration: | 1.8mg/ml | | | | | |
| Buffer: | 10mM PB, pH 7.4 ± 0.2 containing 150mM NaCl | | | | | |
| Preservative: | None | | | | | |
| Applications: | Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00101</td><td>L1G00102</td></tr></table> | | <u>Capture</u> | <u>Detection</u> | L1G00101 | L1G00102 |
| <u>Capture</u> | <u>Detection</u> | | | | | |
| L1G00101 | L1G00102 | | | | | |
| Storage: | Store at <-20°C for three years. | | | | | |