



Hainan Medikamende Trading Co., Ltd.

海南省迈迪卡门德商贸有限责任公司
Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	McAb to Pro-GRP	Catalog #:	L1C03202				
Description:	Mouse Anti Human Pro-GRP Monoclonal Antibody						
Specificity:	Recognizes an epitope on the Human Pro-GRP Antigen.						
Clone:	GR002						
Host Animal:	Mouse. Hybridization of X63-Ag8 myeloma cells with spleen cells from Balb/c mice						
Isotype:	IgG						
Source:	Mouse Ascites						
Immunogen:	Human Pro-GRP Antigen						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration&Lot:	4.2 mg/ml	Lot:	L180201				
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA, CLIA, RIA etc. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C03201</td><td>L1C03202</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C03201	L1C03202
<u>Capture</u>	<u>Detection</u>						
L1C03201	L1C03202						
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY.

海南省迈迪卡门德商贸有限责任公司
地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
邮编：572000
电话：+86-898-31320087，传真：+86-898-31320087
网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.
Address: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,
Xiuyingqu, Haikou, 572000, Hainan, P.R. China
Phone: +86-898-31320087 Fax: +86-898-31320087
Website: www.medical-ivd.com Email: info@medical-ivd.com