

**SPECIFICATION SHEET**

Product Name: McAb to Human SCCA **Catalog #:** L1C01401

Description: Monoclonal Antibody to Squamous cell carcinoma antigen,

Specificity: Recognizes an epitope on the Squamous cell carcinoma antigen.

Clone: SC001

Host Animal: Mouse

Isotype: IgG

Source: Ascites

Immunogen: Squamous cell carcinoma antigen

Format: Purified, Liquid

Purification: Protein A chromatography

Concentration&Lot: 3.3mg/ml **Lot:** L180501

Buffer: 0.01M PBS, 0.15M NaCl, pH 7.4

Preservative: None

Applications: Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Recommended pair for immunoassay:

<u>Capture</u>	<u>Detection</u>
L1C01401	L1C01402

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.