

SPECIFICATION SHEET

Product Name:	McAb to Human SCCA	Catalog #:	L1C01402				
Description:	Monoclonal Antibody to Squamous cell carcinoma antigen,						
Specificity:	Recognizes an epitope on the Squamous cell carcinoma antigen.						
Clone:	SC002						
Host Animal:	Mouse						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Squamous cell carcinoma antigen						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration&Lot:	3.6 mg/ml	Lot:	L180501				
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C01401</td><td>L1C01402</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C01401	L1C01402
<u>Capture</u>	<u>Detection</u>						
L1C01401	L1C01402						

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.