

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Description:</b>	Mouse anti- transferrin, TRF antibody <b>Catalog#:</b> L1K01001						
<b>Background:</b>	Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron (Fe) in biological fluids. Human transferrin is encoded by the TF gene.						
<b>Specificity:</b>	Human						
<b>Clone:</b>	TRF001						
<b>Host Animal:</b>	Mouse.						
<b>Isotype:</b>	IgG						
<b>Immunogen:</b>	Purified Human Serum Albumin						
<b>Format:</b>	Purified, Liquid						
<b>Purification:</b>	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4.						
<b>Concentration&amp;Lot.:</b>	2.4 mg/ml <b>Lot:</b> L190501						
<b>Buffer:</b>	PBS, pH 7.4						
<b>Applications:</b>	Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. The following pair is recommended for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L3I00701</td><td>L1K00401</td></tr><tr><td>L1K00401</td><td>L1K00402, L1K00403</td></tr></table>	<u>Capture</u>	<u>Detection</u>	L3I00701	L1K00401	L1K00401	L1K00402, L1K00403
<u>Capture</u>	<u>Detection</u>						
L3I00701	L1K00401						
L1K00401	L1K00402, L1K00403						
<b>Storage&amp;Shelflife:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw						