

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Description:</b>	Mouse anti- transferrin, TRF antibody				
<b>Background:</b>	Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron (Fe) in biological fluids. Human transferrin is encoded by the TF gene.				
<b>Specificity:</b>	Human				
<b>Clone:</b>	TRF002				
<b>Host Animal:</b>	Mouse.				
<b>Isotype:</b>	IgG				
<b>Immunogen:</b>	Purified Human Serum Albumin				
<b>Format:</b>	Purified, Liquid				
<b>Purification:</b>	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4.				
<b>Concentration&amp;Lot.:</b>	4.2 mg/ml <b>Lot:</b> L181201				
<b>Buffer:</b>	PBS, pH 7.4				
<b>Applications:</b>	Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. The following pair is recommended for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1K01001</td><td>L1K01002</td></tr></table>	<u>Capture</u>	<u>Detection</u>	L1K01001	L1K01002
<u>Capture</u>	<u>Detection</u>				
L1K01001	L1K01002				
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.				