

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

### SPECIFICATION SHEET

**ProductName:** G-17-Ag-BSA **Catalog #:** L5B01001

**Description:** BSA Conjugated Human GASTRIN 17 (G-17)

**Format:** liquid, G-17 purified by HPLC

**Purification:** Purified

**Concentration& Lot.:** 10mg/ml **Lot: L190901**

**Buffer:** 10Mm PBS, pH 7.5

**Preservative:** None

#### Description:

Endogenous peptide, mainly found in the gut, that activates gastric acid secretion, delays gastric emptying and increases blood flow in the gastric mucosa.

**Storage&Shelflife:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.