

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L2K00201	Lot #:	L180501
Description:	BSA-S-adenosyl-L-homocysteine BSA conjugated SAH		
Type:	Synthetic		
Format:	Purified, Liquid		
Concentration:	7.7mg/ml		
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	ELISA		
Storage:	Store at -20°C for long term storage. Avoid repeated freeze-thaw cycles.		