

SPECIFICATION SHEET

Description:	Human CA125 Antigen	Catalog #:	L2C011
Description:	Human Cancer Antigen 125 (CA125) (Ovarian Cancer) (Gastrointestinal Cancer)		
Source:	Human pleural fluids		
Format:	Calibrator grade, Liquid		
Purification:	Gel filtration and column chromatography.		
Concentration:	300KU/ml, (Roche ECLIA)	Lot #:	L160801
Buffer:	20Mm PBS, pH 7.4 ± 0.02		
Preservative:	0.1% Sodium Azide		
Applications:	Specific methodologies have not been tested using this product.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Inactivation:	Not applicable.		
Warnings:	All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		

This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.