



SPECIFICATION SHEET

Product Name: Human CA24-2 Calib. Grade **Catalog #:** L2C00602

Description: Human Cancer Antigen 24-2 (CA24-2) (Gastrointestinal Cancer)
Calibration Grade

Source: Human colion carcinoma.

Format: Purified, Liquid

Purification: Gel filtration and column chromatography.

Concentration: 450KU/ml (Canag ELISA) **Lot #:** L170201

Buffer: PBS, pH 7.4

Preservative: 0.05% Sodium Azide

Applications: Immunogen, calibrator in immunoassays and in vitro studies. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Upon receipt, store at -20°C.

Inactivation: Not applicable.

Warnings: All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.