



SPECIFICATION SHEET

Product Name: Recombinant NSE protein **Catalog #:** L2C00302

Description: Recombinant Neuron Specific Enolase (NSE), Calibrator Grade

Molecular Weight: 56 466 (calculated from sequence)
56 000 (determined by SDS gelelectrophoresis)

Source: E. coli

Format: Purified

Purification: >90% (determined by SDS gelelectrophoresis).

Concentration: Specific protein concentration:1.0mg/ml

Lot #: L190501

Buffer: 10mM Tris-HCL, PH8.0 ± 0.2

Preservative: 0.05% NaN₃

Storage: Store at -20°C. Avoid multiple freeze/thaw cycles.

Warnings:

This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.