

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName: CMV Angiten **Catalog#:** L1Y00201

Description: Cytomegalovirus Ag (P52+PP150)
Cytomegalovirus M Antigen. (P52+PP150)
WM: 25KD

Source: Recombinant ,E.col

Format: Purified, Protein solution in Tris-HCl pH8.0,

Purification: Purity >85%

Concentration&Lot.: 2mg/ml **Lot.:** L180701

Buffer: PBS, pH 7.4 + 0.2 , 0.5mM EDTA

Preservative: None

Applications: Immune colloidal gold immunochromatography

Storage&Shelf life:

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

Warnings: For research only

网址: www.medical-ivd.com 邮箱: info@medical-ivd.com

Website: www.medical-ivd.com Emai: linfo@medical-ivd.com