

海南省迈迪卡门德商贸有限责任公司  
Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

**SPECIFICATION SHEET**

**Catalog #:** L2G00201 **Lot #: L170501**

**Description:** Follicle Stimulating Hormone (FSH) intact human

**Purity:** >85% by SDS-PAGE

**Source:** Human pituitary glands

**Concentration:** 1mg/ml (1800U/mg)

**Contaminants:** not detectable

**Buffer:** 10mM PBS, 150mM NaCl, PH 7.4

**Biohazard:** Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

**Storage:** -20°C or lower, for long term storage.