

---

## SPECIFICATION SHEET

<b>Product Name:</b>	free HCG Beta Antigen	<b>Catalog #:</b>	L2C00501
<b>Description:</b>	Human free Chorionic Gonadotropin beta (free hCG-b)		
<b>Source:</b>	Human Pregnancy Urine		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	≥90% pure (SDS-PAGE).		
<b>Concentration:</b>	1.0 mg/ml	<b>Lot #:</b>	L180301
<b>Buffer:</b>	10mM PBS PH7.4		
<b>Preservative:</b>	0.05%NaN3		
<b>Application:</b>	Specific methodologies have not been tested using this product.		
<b>Storage:</b>	Store at -20°C. Avoid multiple freeze/thaw cycles.		
<b>Inactivation:</b>	Not applicable.		

**Warning:** All human source materials have tested negative for HIV 1, HIV 2, HCV, Hepatitis A antibodies and HBsAg and HIV Antigen. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.