



SPECIFICATION SHEET

Product Name:	Goat anti mouse IgG (H&L)	Catalog #:	L1I00501
Description:	Goat Antibody to mouse Immunoglobulin G (IgG), Heavy & Light chains		
Specificity:	Reacts with heavy chains on mouse IgG and with the light chains on all mouse immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal mouse Serum.		
Host Animal:	Goat		
Isotype:	IgG		
Source:	Serum		
Immunogen:	Purified Mouse IgG		
Purification:	DEAE-chromatography;		
Concentration:	5mg/ml	Lot. :	L180801
Buffer:	0.01M sodium phosphate buffer, 0.15M NaCl, pH 7.2.		
Preservative:	None		
Applications:	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage & Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY