

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Catalog #:</b>	L1G00302	<b>Lot #:</b>	L180801
<b>Description:</b>	PcAb to Thyroxine (T4) Goat anti Thyroxine (T4) Ponoclonal Antibody		
<b>Specificity:</b>	<u>Cross-reactivity:</u> T3<1%		
<b>Clone:</b>	N/A		
<b>Host Animal:</b>	Goat		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	T4-BSA		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure. Chromatographic methods		
<b>Concentration:</b>	1.7mg/ml		
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Specific methodologies have not been tested using this product.		
<b>Storage:</b>	store at <-20°C.		