

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00202	Lot #:	L180801
Description:	PAb to Triiodothyronine (T3) Goat anti Triiodothyronine (T3) Polyclonal Antibody		
Specificity:	<u>Cross-reactivity:</u> T4 <1%		
Clone:	N/A		
Host Animal:	Goat		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	T3-BSA		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE).		
Concentration:	1.6mg/ml		
Buffer:	10mM Phosphate, pH 7.4 containing 150mM Sodium chloride		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage:	Store at <-20°C for at least three years.		