



## SPECIFICATION SHEET

<b>Product Name:</b>	HBc Ag	<b>Catalog #:</b>	L1H00701
<b>Description:</b>	Recombinant Hepatitis B c Ag Hepatitis B c Antigen (Hepatitis B core Antigen, 144aa,14kD)		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure (SDS-PAGE) Ion exchange chromatography		
<b>Concentration:</b>	2.6 mg/ml	<b>Lot.:</b>	L160201
<b>Buffer:</b>	0.02M Sodium phosphate, pH 7.2		
<b>Preservative:</b>	none		
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

### For Research Use Only