



Specification Sheet

Product Name:	HCV MP Antigen	Catalog # :	L1H00801
Description:	Hepatitis C Virus HCV (Core+NS3+NS4+NS5) Antigen, genotypecore+3+4+5.		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	95% pure (SDS-PAGE)		
Concentration:	1.6mg/ml (OD280nm).(Lot. specific)		
Buffer:	10Mm PBS PH7.4		
Preservative:	none		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	All human source materials have tested negative for HIV1, HIV2 and HCV antibodies. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

For Research Use Only.