

SPECIFICATION SHEET

Product Name: HIV-1 P24 **Catalog #:** L2H002

Description: HIV-1 P24 Recombinant antigen, Antigen Grade

Source: Recombinant ,E.coli

Format: Purified, Protein solution in Tris-HCl pH8.0,

Purification: Purity >90%

Concentration&Lot: 2mg/ml **Lot:** L181201

Buffer: 20Mm Tris-HCl, pH 8.0 ± 0.02,

Preservative: None

Application: Immunoassays or Immunization.

Storage&Shelf life: Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.

Warnings: For research only