



SPECIFICATION SHEET

Product Name:	Recombinant HP UreaseA	Catalog #:	L1H01203
Description:	Rcombinant Helicobacter pylori (Hp) Urease A Antigen, Antigen Grade		
Molecular Weight:	26.71KD (calculated from sequence) 26.7KD (determined by SDS gelelectrophoresis)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	> 85% (determined by SDS gelelectrophoresis)		
Concentration&Lot:	1.29mg/ml,	Lot:	L180401
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.1		
Preservative:	0.2%NaN3		
Application:	Immunoassays or Immunization.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		