

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Catalog #:** L2G01201                      **Lot #: L160601**

**Description:** Human Myeloperoxidase purified (MPO)

**Grade:** Affinity purified

**Purity:** 90% by SDS-PAGE

**Source:** Human leukocyte mass

**Concentration:** 3.3mg/ml

**Form:** 20mM Tris-HCl, pH 7.5, 0.5M NaCl

**Biohazard:** Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

**Storage:** -20°C or lower, for long term storage.