

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #: L2G00801 **Lot #:** L161001

Description: Human Thyrogloblin (TG) Purified

Grade: High Purity

Purity: >80% by SDS-PAGE; Ion-exchange chromatography

Source: Thyroid tissue, surgery

Concentration: 6.2mg/ml

Buffer: 10mM PBS, pH 7.4, 0.05% NaN₃.

Biohazard: Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

Storage: Store at -20°C.