

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	Trophoblast Antigen	Catalog #:	L1F004
Description:	Human Trophoblast Antigen (TA-Ag)		
Source:	Human Trophoblast		
Format:	Ion exchange, chromatography, Purified, Solubility is limited at neutral pH and physiological salt concentration.		
Concentration& Lot.:	1.41mg/ml(Commassie Blue Staining Kit)	Lot:	L160701
Buffer:	20mM Tris-Hcl 150Mm NaCl ,PH7.5		
Preservative:	0.2% NaN3		
Applications:	Suitable for use as a standard in immunoassay, Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	A blood sample form the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		