

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name: Zona Pellucid **Catalog #:** L1F005

Description: Human zona pellucid Antigen (ZP-Ag)

Source: Human zona pellucid

Format: Solubility is limited at neutral pH and physiological salt concentration.

Concentration& Lot.: 2.9 mg/ml (Commassie Blue Staining Kit)

Lot: L160301

Buffer: 20mM Tris-Hcl 150Mm NaCl ,PH7.5

Preservative: 0.2% NaN₃

Applications: Suitable for use as a standard in immunoassay,
Each laboratory should determine an optimum working titer for use in its particular application.

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

Warnings: A blood sample form the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.