

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Product Name:** Recombinant IGF2

**Catalog #:** L4I00702

#### Intorduction

Insulin like growth factor 2 (IGF2 )gene was the earliest discovered endogenous imprinted gene, its expression, regulation and control does not follow the Mendel rule. IGF2 is a multifunctional regulatory factor that plays an important role in cell differentiation, proliferation, embryonic growth and tumor cell proliferation.

**Molecular weight:** 15KD

**Source:** E. coli

**Format:** Purified, Liquid

**Purification:** > 90% ,Sephrose derived purification

**Concentration& Lot.:** 1.86 mg/ml      **Lot:** L181202

**Buffer:** 10mm PBS,0.1%SKL

**Preservative:** 0.1%SKL

**Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw