

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Description: L-T₄ L-Thyroxine Sodium Salt Pentahydrate

Catalog #: L2G01002

Molecular Formula: C₁₅H₁₀I₄NO₄Na·5H₂O

MW: 888.93

Purity: 80%

Lot: L160401

Forma: Powder

Storage: Store at 2 ~ 8 °C, Against light or wet.

Dissolved:

Dissolve one by one drops of NaOH 0.1mol/L, and make it dissolve completely.
The solution is diluted to the appropriate concentration with 50% propylene glycol solution. (Note: this method is only applicable to the calibration of the product, for reference. Other uses of the dissolution method please refer to the relevant information.)