

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Catalog #:** L2G00301

**Lot #: L180501**

**Description:** Human Luteinizing Hormone (LH)

**Grade:** Standard

**Purity:** About 85% by SDS-PAGE

**Source:** Human pituitary glands

**Concentration:** 1mg/ml (1800U/mg)

**Form:** 10mM PBS, pH 7.4

**Biohazard:** Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

**Storage:** -20°C or lower, for long term storage.