

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Product Name:</b>	McAb to FFN	<b>Catalog #:</b>	L1F00101
<b>Description:</b>	Monoclonal Antibody to Fetal Fibronectin protein, FFN		
<b>Clone:</b>	FN001		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Purified Fetal Fibronectin protein, FFN		
<b>Format:</b>	Purified, Liquid		
<b>Concentration&amp;Lot:</b>	6.8 mg/ml (OD280nm) <b>Lot #:</b> L180601		
<b>Buffer:</b>	10mM PBS pH 7.4 ± 0.2		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay:		
<b>Storage:</b>	Store at -20°C for at least 3 years.		