



---

## SPECIFICATION SHEET

<b>Catalog #:</b>	L1K00103	<b>Lot #:</b>	L180301
<b>Description:</b>	McAb to Cystatin C antibody Mouse Anti Human Cystatin C Monoclonal Antibody		
<b>Specificity:</b>	Recognizes an epitope on Cystatin C Antigen		
<b>Clone:</b>	CY003		
<b>Host Animal:</b>	Mouse		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Recombinant Cystatin C Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	Protein A chromatography		
<b>Concentration:</b>	3.8mg/ml		
<b>Buffer:</b>	0.01M PBS, 0.15M NaCl, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in WB, IHC		
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		