

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	McAb to Human Albumin	Catalog #:	L3I00403
Description:	Mouse anti- Human Serum Albumin (HSA) antibody		
Background:	Human serum albumin is the most abundant protein in human blood plasma. It is produced in the liver. Albumin comprises about half of the blood serum protein. It is soluble and monomeric.		
Specificity:	Human		
Clone:	HA003		
Host Animal:	Mouse.		
Isotype:	IgG		
Immunogen:	Purified Human Serum Albumin		
Format:	Purified, Liquid		
Purification:	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4. No preservatives added.		
Concentration&Lot.:	3.3 mg/ml (Lot. specific)		
Buffer:	PBS, pH 7.4		

海南省迈迪卡门德商贸有限责任公司

地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502

邮编: 572000

电话: +86-898-31320087

传真: +86-898-31320087

网址: www.medical-ivd.com

邮箱: info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.

Adress: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

Phone: +86-898-31320087

Fax: +86-898-31320087

Website: www.medical-ivd.com

Emai: linfo@medical-ivd.com

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Applications: Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. The following pair is recommended for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
L3I00401	L3I00402, L3I00403

Storage& Shelf life: Store at -20°C . Avoid multiple freeze/thaw cycles.

海南省迈迪卡门德商贸有限责任公司

地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502

邮编: 572000

电话: +86-898-31320087 传真: +86-898-31320087

网址: www.medical-ivd.com

邮箱: info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.

Adress: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

Phone: +86-898-31320087 Fax: +86-898-31320087

Website: www.medical-ivd.com

Emai: linfo@medical-ivd.com