



SPECIFICATION SHEET

Product Name:	Mycoplasma pneumoniae MP P1Ag	Catalog #:	L3H002
Description:	Recombinant Mycoplasma pneumoniae (MP-Ag) Antigen Grade		
Source:	Recombinant,E.coli		
Format:	Purified, Protein solution in PBS		
Purification:	Purity >90%		
Concentration&Lot:	2.4 mg/ml,	Lot:	L190601
Buffer:	10Mm PBS, 2Mm EDTA, pH 7.4± 0.02		
Preservative:	None		
Applications:	Immune colloidal gold immunochromatography		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		