



SPECIFICATION SHEET

Catalog #:	L2L009	Lot #:	L190601
Description:	Procollagen Type III protein		
Source:	Human placenta		
Background:	Collagens are the major fibrillar components of most connective tissues. Procollagens are Interstitial collagens synthesised and secreted out of cells. Type III procollagen (PIIINP) is a sensitive marker of liver and lung fibrosis.		
Purification:	purified.		
Concentration:	2.36mg/ml		
Buffer:	0.5M CH ₃ COOH		
Preservative:	None		
Applications:	User optimized		
Storage:	Store at -20 °C. Avoid multiple freeze/thaw cycles.		