



## SPECIFICATION SHEET

<b>Product Name:</b>	Rabbit anti FITC	<b>Catalog #:</b>	L1101402
<b>Description:</b>	Rabbit Antibody to fluorescein isothiocyanate, FITC PcAb		
<b>Specificity:</b>	Reacts with fluorescein isothiocyanate, FITC		
<b>Host Animal:</b>	Rabbit		
<b>Immunogen:</b>	Purified BSA conjugated FITC		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure (SDS-PAGE). Purified.		
<b>Concentration:</b>	1.2mg/ml	<b>Lot:</b>	L180601
<b>Buffer:</b>	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2		
<b>Preservative:</b>	None		
<b>Applications:</b>	Specific methodologies have not been tested using this product.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

**FOR RESEARCH USE ONLY**