

SPECIFICATION SHEET

Product Name: Sheep IgG **Catalog #:** L11019

Description: Sheep ImmunoglobulinG(IgG)

Isotype: IgG

Source: Normal Sheep Serum..

Immunogen: Purified Sheep IgG

Format: 0.01M PBS pH 7.4

Concentration: 25mg/ml **Lot:** L171001

Purification: Octoic acid purified.

Application: Sepecific methodologies have not been tested using this product.

Preservative: None

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY