



---

## SPECIFICATION SHEET

**Product Name:** Rabbit anti Chicken IgY **Catalog #:** L1102001

**Description:** Rabbit anti Chicken Immunoglobulin Y Polyclonal antibody

**Specificity:** Reacts with heavy chains on Immunoglobulin Y (abbreviated as IgY).

**Isotype:** IgG

**Source:** Rabbit

**Format:** 0.01M PBS with 0.15M NaCl pH 7.4

**Concentration:** 29mg/ml **Lot:** L190301

**Purification:** Purified

**Application:** User optimized

**Preservative:** None

**Storage & Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

**FOR RESEARCH USE ONLY**