

SPECIFICATION SHEET

Description:	Prostate Specific Antigen (PSA)	Catalog #:	L2C001
Source:	Human Seminal Fluid		
Format:	Purified, Liquid		
Activity Concentration&Lot:	1 mg/ml (Canag ELISA)	Lot:	L181101
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4		
Preservative:	none		
Application:	Specific methodologies have not been tested using this product.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		

All donor units tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.

FOR RESEARCH USE ONLY.