

SPECIFICATION SHEET

Product Name:	Purified D-dimmer Antigen	Catalog #:	L4C00801
Description:	Purified native Human D-Dimer protein		
Source:	Human serum		
Purification:	Affinity purified.		
Concentration&Lot:	5mg/ml	Lot:	L20200501
Buffer:	pH 7.5 PBS 150Mm NaCl		
Preservative:	None		
Applications:	Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.