

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L2G00901	Lot #:	L180901
Description:	Purified Fetal Fibronectin protein		
Background:	FF N, a molecular weight of 310 ~ 550 KD, composed of two polypeptide chains, containing a specific epitope (m CS zone) in which the fetus tumor area.		
Source:	Human placenta		
Format:	Purified, Liquid		
Purification:	About $\geq 90\%$ pure. Affinity purified.		
Concentration:	2 mg/ml		
Buffer:	20mmNH ₄ AC PH5.0		
Storage:	Aliquot and store at -20°C. Avoid repeated Freeze/Thaw cycles		

Biohazard Information:

Donor samples were tested and found to be negative for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II. Nonetheless caution should be used when handling this material as there is a margin of error in all tests.