



SPECIFICATION SHEET

Catalog #:	L4I01102	Lot #:	L180701
Description:	Purified RBP-4 (Retinol-Binding Protein 4,RBP 4)		
Source:	Human, Native		
Format:	Purified, Lyophilized		
Purification:	>90% (determined by SDS gelelectrophoresis).		
Concentration:	Lypholized,Reconstitute by adding 1ml DD water to 1.0mg/ml		
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Storage:	Store at -20°C. Avoid multiple freeze/thaw cycles.		
Warnings:	All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.